

WEAVER
and company

Nūprep®

Abrasive Skin Prep Gel



*Avoid Eye Contact

Nuprep[®] Skin Prepping Gel

Nuprep is a mild abrasive gel that strips away dry surface skin and moistens the corneum skin layer. Gently applied, Nuprep Gel improves conductivity and helps achieve high quality tracings.

Nuprep Gel has earned an excellent reputation through extensive use by leading cardio-technologists in medical facilities throughout the world.

Nuprep Gel is valuable whenever a reduction in skin impedance would enhance performance of a monitoring electrode. Gently apply Nuprep Gel to the skin site prior to attaching the electrode. Nuprep Gel virtually eliminates problems such as diaphoresis and muscle artifact.

Zero in on these procedures with Nuprep Gel.

- Cardiac Stress Test
- Holter Monitoring
- Cardiac Rehabilitation Monitoring
- Cardiac Catheter Monitoring
- ECG Telemetry



4oz (114 Gm) Tube, 3-pack
Product # 10-30



25 Gm Tube, 6-pack
Product # 10-61

FOR MORE INFORMATION
Contact:

Weaver and Company
565 Nucla Way Unit B
Aurora, CO 80011
Toll Free 800-525-2130
Fax 303-367-5118
www.doweaver.com

WEAVER
and company